

# SimpleChIP® Enzymatic Chromatin IP Kit (Agarose Beads)

✓ 1 Kit  
(30 Immunoprecipitations)



**Orders** ■ 877-616-CELL (2355)  
orders@cellsignaling.com  
**Support** ■ 877-678-TECH (8324)  
info@cellsignaling.com  
**Web** ■ www.cellsignaling.com

rev. 05/31/12

**For Research Use Only. Not For Use In Diagnostic Procedures.**

Box 1 contains the following buffers/reagents:	
Components to be stored at 4°C:	
Glycine Solution (10X)	100 ml
Buffer A (4X)	25 ml
Buffer B (4X)	25 ml
ChIP Buffer (10X)	20 ml
ChIP Elution Buffer (2X)	7 ml
5 M NaCl	3 ml
0.5 M EDTA	1 ml
ChIP-Grade Protein G Agarose Beads (blocked with BSA and sonicated salmon sperm DNA)	1 ml
Components to be stored at room temperature:	
DNA Binding Reagent A (add 12 ml isopropanol before use)	12 ml
DNA Wash Reagent B (add 24 ml ethanol before use)	6 ml
DNA Elution Reagent C	2 X 1 ml
DNA Spin Columns	36

Box 2 contains the following buffers/reagents:	
Components to be stored at -20°C:	
Protease Inhibitor Cocktail (200X)	500 µl
RNAse A (10 mg/ml)	50 µl
Micrococcal Nuclease (2000 gel units/µl)	60 µl
Proteinase K (20 mg/ml)	100 µl
SimpleChIP® Human RPL30 Exon 3 Primers 1	150 µl
SimpleChIP® Mouse RPL30 Intron 2 Primers 1	150 µl
Histone H3 (D2B12) XP® Rabbit mAb (ChIP Formulated) #4620	100 µl
Normal Rabbit IgG (1 µg/µl) #2729	50 µl
1M DTT	200 µl

**Storage:** Please store components at the temperatures indicated on the individual tube labels.

**Note:** Buffer A (4X), Buffer B (4X) and ChIP-Grade Protein G Agarose Beads contain 0.05% sodium azide.


## Reagents not Supplied:

PMSF (0.1 M stock)  
Formaldehyde (37% Stock)  
Ethanol (96-100%)  
Isopropanol  
1X PBS  
Nuclease-free water  
Taq DNA polymerase  
dNTP mix


Please visit [www.cellsignaling.com/technologies/chip.html](http://www.cellsignaling.com/technologies/chip.html) for a complete listing of recommended companion products.

**Description:** The SimpleChIP® Enzymatic Chromatin IP Kit provides all the reagents required for performing 6 chromatin preparations (or optimizations) and 30 chromatin immunoprecipitation assays and is optimized for 4 X 10<sup>7</sup> cells per experiment. Kit components are stable for 1 year from date of shipment when stored as directed.


**Specificity/Sensitivity:** The SimpleChIP® Enzymatic Chromatin IP Kit can be utilized with any ChIP-validated antibody to detect endogenous levels of protein-DNA interactions and histone modifications in mammalian cells (see Figures 2 and 3). The positive control Histone H3 Rabbit mAb recognizes many different species of the highly conserved Histone H3 protein, including human, mouse, rat and monkey. Primer sets are included for the human and mouse positive control RPL30 gene loci; however, the use of other species with the kit requires the design of additional control primer sets.



**FIGURE 2.** Chromatin immunoprecipitations were performed using digested chromatin from HeLa cells and either Histone H3 (D2B12) XP® Rabbit mAb (ChIP Formulated) #4620 (lane 2), Rpb1 CTD (4H8) Mouse mAb #2629 (lane 3), Di-Methyl Histone H3 (Lys9) Antibody #9753 (lane 4), or Normal Rabbit IgG #2729 (lane 5). Purified DNA was analyzed by standard PCR methods using SimpleChIP® Human RPL30 Exon 3 Primers #7014, SimpleChIP® Human MyoD1 Exon 1 Primers #4490, and SimpleChIP® Human  $\alpha$  Satellite Repeat Primers #4486. PCR products were observed for each primer set in the input sample (lane 1) and various ChIP samples, but not in the Normal Rabbit IgG ChIP sample (lane 5).



**FIGURE 3.** Chromatin immunoprecipitations were performed using digested chromatin from HeLa cells and the indicated ChIP-validated antibodies. Purified DNA was analyzed by quantitative Real-Time PCR, using SimpleChIP® Human RPL30 Exon 3 Primers #7014 (control primer set), SimpleChIP® Human MyoD1 Exon 1 Primers #4490, and SimpleChIP® Human  $\alpha$  Satellite Repeat Primers #4486. The amount of immunoprecipitated DNA in each sample is represented as signal relative to the total amount of input chromatin (equivalent to 1).



**FIGURE 1:** HeLa cells were formaldehyde-crosslinked and chromatin was prepared and digested as described in Section A of protocol. DNA was purified as described in Section B and 10 µl were separated by electrophoresis on a 1% agarose gel (lane 2) and stained with ethidium bromide. Lane 2 shows that the majority of chromatin was digested to 1 to 5 nucleosomes in length (150 to 900 bp).

SimpleChIP® kit co-developed by CST and NEB; Antibodies provided by CST; Enzymes provided by NEB.



## Introduction:

The chromatin immunoprecipitation (ChIP) assay is a powerful and versatile technique used for probing protein-DNA interactions within the natural chromatin context of the cell (1,2). This assay can be used to identify multiple proteins associated with a specific region of the genome, or the opposite, to identify the many regions of the genome associated with a particular protein (3-6). In addition, the ChIP assay can be used to define the spatial and temporal relationship of a particular protein-DNA interaction. For example, the ChIP assay can be used to determine the specific order of recruitment of various protein factors to a gene promoter or to "measure" the relative amount of a particular histone modification across an entire gene locus during gene activation (3,4). In addition to histone proteins, the ChIP assay can also be used to analyze binding of transcription factors, transcription co-factors, DNA replication factors and DNA repair proteins.

When performing the ChIP assay, cells are first fixed with formaldehyde, a reversible protein-DNA cross-linking agent that serves to fix or "preserve" the protein-DNA interactions occurring in the cell (see method overview) (1,2). Cells are then lysed and chromatin is harvested and fragmented using either sonication or enzymatic digestion. The chromatin is then subjected to immunoprecipitation using antibodies specific to a particular protein or histone modification. Any DNA sequences that are associated with the protein or histone modification of interest will co-precipitate as part of the cross-linked chromatin complex and the relative amount of that DNA sequence will be enriched by the immunoselection process. After immunoprecipitation, the protein-DNA cross-links are reversed and the DNA is purified. The enrichment of a particular DNA sequence or sequences can then be detected by a number of different methods.

Standard PCR methods are often employed to identify the DNA sequences or regions of the genome associated with a particular protein or histone modification (1,2). PCR is used to measure the relative abundance of a particular DNA sequence enriched by a protein-specific immunoprecipitation versus an immunoprecipitation with a non-specific antibody control. PCR products are run on an agarose or acrylamide gel to facilitate

quantification, and the level of enrichment of the DNA sequence is determined relative to the total amount of input DNA (percent of input). The level of enrichment can also be expressed as fold enrichment above background (enrichment relative to that of the non-specific antibody control). Real-Time PCR provides a more accurate, gel-free system for the quantification of DNA enrichment. Alternatively, the ChIP assay can be combined with genomic tiling micro-array (ChIP on chip) techniques, sequencing, or cloning strategies, which allow for genome-wide analysis of protein-DNA interactions and histone modifications (5-8).

The SimpleChIP® Kit contains buffers and reagents needed to perform the ChIP assay with mammalian cells and works for both histone modifications and non-histone DNA-binding proteins. After cell lysis, the chromatin is fragmented by partial digestion with Micrococcal Nuclease to obtain chromatin fragments of 1 to 5 nucleosomes in size. Enzymatic digestion of chromatin is much milder than sonication and eliminates problems due to variability in sonication power and emulsification of chromatin during sonication, which can result in incomplete fragmentation of chromatin or loss of antibody epitopes due to protein denaturation and degradation. The chromatin immunoprecipitations are performed using antibodies and either ChIP Grade Protein G Agarose or ChIP Grade Protein G Magnetic Beads. After reversal of protein-DNA cross-links, the DNA is purified using DNA purification spin columns provided in the kit. The DNA purification spin columns combine the convenience of spin-column technology with the selective binding properties of a uniquely designed silica membrane that allows for efficient recovery of DNA and removal of protein contaminants without the need for phenol/chloroform extractions and ethanol precipitations. After DNA purification, the enrichment of particular DNA sequences can be analyzed by a variety of methods, including standard PCR, quantitative real-time PCR, amplification for ChIP on chip, sequencing or cloning techniques.

In addition to providing buffers and reagents required to perform the ChIP assay, the SimpleChIP® Kit provides important controls that allow for user determination of a successful ChIP experiment. The kit contains a positive control Histone H3

Rabbit mAb, a negative control Normal Rabbit IgG Antibody and primer sets for PCR detection of the ribosomal protein L30 (RPL30) gene locus (human and mouse primer sets included). Histone H3 is a core component of chromatin in the cell and is bound to most DNA sequences throughout the genome, including the RPL30 locus. Thus, immunoprecipitation of chromatin with the Histone H3 antibody will enrich for the RPL30 gene, while immunoprecipitation with the Normal Rabbit IgG will not result in RPL30 gene enrichment. This enrichment can be quantified using either standard PCR or quantitative real-time PCR methods and the RPL30 primer sets provided in the kit. Importantly, since histone H3 is bound to most DNA sequences throughout the genome, the Histone H3 Rabbit mAb serves as a positive control IP for almost any locus studied, giving the user even more confidence that their ChIP experiment was performed successfully.

The SimpleChIP® Kit provides enough reagents to perform up to 6 chromatin preparations (or optimizations) and 30 immunoprecipitations and is optimized for  $4 \times 10^7$  cells per experiment. A ChIP assay can be performed in as little as two days and can easily be scaled up or down for use with more or fewer cells.

## Background References:

- (1) Orlando, V. (2000) *Trends Biochem Sci* 25, 99-104.
- (2) Kuo, M.H. and Allis, C.D. (1999) *Methods* 19, 425-33.
- (3) Agalioti, T. et al. (2000) *Cell* 103, 667-78.
- (4) Soutoglou, E. and Talianidis, I. (2002) *Science* 295, 1901-4.
- (5) Mikkelsen, T.S. et al. (2007) *Nature* 448, 553-60.
- (6) Lee, T.I. et al. (2006) *Cell* 125, 301-13.
- (7) Weinmann, A.S. and Farnham, P.J. (2002) *Methods* 26, 37-47.
- (8) Wells, J. and Farnham, P.J. (2002) *Methods* 26, 48-56.


## Method Overview

Cells are fixed with formaldehyde to cross-link histone and non-histone proteins to DNA.

Chromatin is digested with Micrococcal Nuclease into 150-900 bp DNA/protein fragments.

Antibodies specific to histone or non-histone proteins are added and the complex co-precipitates and is captured by Protein G Agarose or Protein G magnetic beads.

Cross-links are reversed, and DNA is purified and ready for analysis.



# Chromatin Immunoprecipitation Protocol

## A *In Vivo* Crosslinking, Nuclei Preparation and Nuclease S7 Digestion of Chromatin:

### Before starting:

Stimulate or treat approximately  $4 \times 10^7$  cells for each experiment. This number of cells will generate one chromatin preparation that can be used for up to 10 separate immunoprecipitations.

- For HeLa cells, this is equivalent to 5 x 15 cm culture dishes containing cells that are 90% confluent in 20 ml of growth media.
- Include one extra dish to be used for determination of cell number using a hemocytometer.

### Remove and warm:

- 200X Protease Inhibitor Cocktail
- 1 M DTT
- 0.1 M PMSF (not supplied)
- 10X ChIP Buffer and ensure SDS is completely in solution.

### Prepare and chill on ice for each set of $4 \times 10^7$ cells (5 dishes):

- 10 ml 10X glycine
- 200 ml 1X PBS (not supplied)
- 10 ml 1X PBS + 100  $\mu$ l PMSF
- 10 ml 1X Buffer A (2.5 ml 4X Buffer A + 7.5 ml water) + 5  $\mu$ l 1M DTT + 50  $\mu$ l 200X Protease Inhibitor Cocktail (PIC) + 100  $\mu$ l PMSF
- 11 ml 1X Buffer B (2.75 ml 4X Buffer B + 8.25 mL water) + 5.5  $\mu$ l 1 M DTT
- 1 ml 1X ChIP Buffer (100  $\mu$ l 10X ChIP Buffer + 900  $\mu$ l water) + 5  $\mu$ l 200X Protease Inhibitor Cocktail (PIC) + 10  $\mu$ l PMSF.

1. To crosslink proteins to DNA, add 540  $\mu$ l of 37% formaldehyde to each 15 cm culture dish containing 20 ml medium. Swirl briefly to mix and incubate for 10 minutes at room temperature.
  - Final formaldehyde concentration is 1%.
  - Use fresh formaldehyde that is not past the manufacturer's expiration date.
  - Addition of formaldehyde may result in a color change of the medium
2. Add 2 ml of 10X glycine to each 15 cm dish, swirl briefly to mix, and incubate for 5 minutes at room temperature.
  - Addition of glycine may result in a color change of the medium.
3. For suspension cells, transfer cells to a 50 ml conical tube, centrifuge at 1,500 rpm in a bench top centrifuge for 5 min at 4°C, and wash pellet two times with 20 ml ice-cold 1X PBS. Immediately proceed to step 5.
4. For adherent cells, remove media and wash cells two times with 20 ml ice-cold 1X PBS, completely removing wash from culture dish each time. Add 2 ml ice-cold 1X PBS + PMSF to the 15 cm dish of cells. Scrape and resuspend cells into cold buffer. Combine cells from all 5 plates into one 15 ml conical tube. Centrifuge cells at 1,500 rpm in a bench top centrifuge for 5 min at 4°C.
5. Remove supernatant and resuspend cells in 10 ml ice-cold Buffer A + DTT + PIC + PMSF. Incubate on ice for 10 min. Mix by inverting tube every 3 min.
6. Pellet nuclei by centrifugation at 3,000 rpm in a bench top centrifuge for 5 min at 4°C. Remove supernatant and resuspend pellet in 10 ml ice-cold Buffer B + DTT. Repeat centrifugation, remove supernatant, and resuspend pellet in 1.0 ml Buffer B + DTT. Transfer sample to a 1.5 ml microcentrifuge tube.
7. Add 5  $\mu$ l of Micrococcal Nuclease, mix by inverting tube several times and incubate for 20 min at 37°C with frequent mixing to digest DNA to length of approximately 150-900 bp. Mix by inversion every 3 to 5 min.
  - The amount of Micrococcal Nuclease required to digest DNA to the optimal length may need to be determined empirically for individual cell lines (see Appendix A).
  - HeLa nuclei digested with 5  $\mu$ l Micrococcal Nuclease per  $4 \times 10^7$  cells in 1 ml Buffer B + DTT gave the appropriate length DNA fragments (see Figure 1).

8. Stop digest by adding 100  $\mu$ l of 0.5 M EDTA and placing tube on ice.
9. Pellet nuclei by centrifugation at 13,000 rpm in a microcentrifuge for 1 min at 4°C and remove supernatant.
10. Resuspend nuclear pellet in 1 ml of 1X ChIP buffer + PIC + PMSF and split into two tubes of 500  $\mu$ l. Incubate on ice for 10 min.
11. Sonicate each tube of lysate with several pulses to break nuclear membrane. Incubate samples for 30 sec on wet ice between pulses.
  - Optimal conditions required for complete lysis of nuclei can be determined by observing nuclei under light microscope before and after sonication.
  - HeLa nuclei were completely lysed after 3 sets of 20-second pulses using a VirTis Virsonic 100 Ultrasonic Homogenizer/Sonicator at setting 6 with a 1/8-inch probe.
  - Alternatively, nuclei can be lysed by homogenizing the lysate 20 times in a Dounce homogenizer; however, lysis may not be as complete.
12. Clarify lysates by centrifugation at 10,000 rpm in a microcentrifuge for 10 min at 4°C.
13. Transfer supernatant to a new tube. This is the cross-linked chromatin preparation, which should be stored at -80°C until further use. Remove 50  $\mu$ l of the chromatin preparation for Analysis of Chromatin Digestion and Concentration (Section B).

## B Analysis of Chromatin Digestion and Concentration (Recommended Step)

1. To the 50  $\mu$ l chromatin sample (from Step 13 in Section A), add 100  $\mu$ l nuclease-free water, 6  $\mu$ l 5 M NaCl and 2  $\mu$ l RNase A. Vortex to mix and incubate samples at 37°C for 30 min.
2. To each RNase A-digested sample, add 2  $\mu$ l Proteinase K. Vortex to mix and incubate samples at 65°C for 2 hours.
3. Purify DNA from samples using spin columns as described in Section F.
4. After purification of DNA, remove a 10  $\mu$ l sample and determine DNA fragment size by electrophoresis on a 1% agarose gel with a 100 bp DNA marker. DNA should be digested to a length of approximately 150-900 bp (1 to 5 nucleosomes; see Figure 1).
5. To determine DNA concentration, transfer 2  $\mu$ l of purified DNA to 98  $\mu$ l TE to give a 50-fold dilution and read the OD<sub>260</sub>. The concentration of DNA in  $\mu$ g/ml is OD<sub>260</sub> x 2,500. DNA concentration should ideally be between 50 and 200  $\mu$ g/ml.

**NOTE:** For optimal ChIP results, it is highly critical that the chromatin is of appropriate size and concentration. Over-digestion of chromatin may diminish signal in the PCR quantification. Under-digestion of chromatin may lead to increased background signal and lower resolution. Adding too little chromatin to the IP may result in diminished signal in the PCR quantification. A protocol for optimization of chromatin digestion can be found in Appendix A.

## C Chromatin Immunoprecipitation:

### Before starting:

### Remove and warm:

- 200X Protease Inhibitor Cocktail
- 10X ChIP Buffer and ensure SDS is completely in solution.

### Remove and place on ice:

- Cross-linked chromatin preparation (from Step 13 in Section A)
- ChIP-Grade Protein G Agarose Beads #9007
- Antibodies for IP. This includes the positive control Histone H3 (D2B12) XP® Rabbit mAb (ChIP Formulated) #4620 and negative control Normal Rabbit IgG #2729.

## Chromatin Immunoprecipitation Protocol (cont.)

**NOTE:** The cross-linked chromatin preparation does not need to be diluted as described below. Antibodies can be added directly to the undiluted chromatin preparation for immunoprecipitation of chromatin complexes.

1. In one tube, prepare enough 1X ChIP Buffer for the desired number of immunoprecipitations. Each precipitation contains 400 µl of 1X ChIP Buffer (40 µl of 10X ChIP Buffer + 360 µl water) and 2 µl Protease Inhibitor Cocktail. When determining the number of immunoprecipitations, the user should include the positive control Histone H3 (D2B12) XP® Rabbit mAb (ChIP Formulated) #4620 and negative control Normal Rabbit IgG #2729 samples. Place mix on ice.
2. To the prepared ChIP buffer, add the equivalent of 100 µl (5 to 10 µg of chromatin DNA) of the cross-linked chromatin preparation (from Step 13 in Section A) for each immunoprecipitation. For example, for 10 immunoprecipitations, prepare a tube containing 4 ml 1X ChIP Buffer (400 µl 10X ChIP Buffer + 3.6 ml water) + 20 µl Protease Inhibitor Cocktail + 1 ml digested chromatin preparation.
3. Remove a 10 µl sample of the diluted chromatin and transfer to a microfuge tube. This is your 2% Input Sample, which can be stored at -20°C until further use (Step 1 in Section E).
4. For each immunoprecipitation, transfer 500 µl of the diluted chromatin to a microfuge tube and add the immunoprecipitating antibody. The amount of antibody required for IP varies for each antibody and should be determined by the user. For the positive control Histone H3 (D2B12) XP® Rabbit mAb (ChIP Formulated) #4620, add 10 µl to the IP sample. For the negative control Normal Rabbit IgG #2729, add 1 µl (1 µg) to 5 µl (5 µg) to the IP sample. Incubate IP samples 4 hours to overnight at 4°C with rotation.
5. Add 30 µl of ChIP-Grade Protein G Agarose Beads #9007 and incubate for 2 h at 4°C with rotation.
6. Proceed to Section D.

### D Washing of the Immunoprecipitated Chromatin:

**Before starting:**

**Remove and warm:**

- 10X ChIP Buffer and ensure SDS is completely in solution.

**Prepare and chill on ice for each IP:**

- **Low salt wash:** 3 ml 1X ChIP Buffer (300 µl 10X ChIP Buffer + 2.7 ml water)
- **High salt wash:** 1 ml 1X ChIP Buffer (100 µl 10X ChIP Buffer + 900 µl water) + 70 µl 5M NaCl

1. Pellet Protein G Agarose Beads in each immunoprecipitation (from Step 5, Section C) by brief 1 min centrifugation at 6000 rpm and remove supernatant.
2. Add 1 ml of low salt wash to the beads and incubate at 4°C for 5 min with rotation. Repeat steps 1 and 2 two additional times, for a total of 3 low salt washes.
3. Add 1 ml of high salt wash to the beads and incubate at 4°C for 5 min with rotation. Immediately proceed to Section E.

### E Elution of Chromatin from Antibody/Protein G Beads and Reversal of Cross-links:

**Before starting:**

- Remove and warm 2X ChIP Elution Buffer in a 37°C water bath and ensure SDS is in solution.
- Set a water bath or thermomixer to 65°C.
- Prepare 150 µl 1X ChIP Elution Buffer (75 µl 2X ChIP Elution Buffer + 75 µl water) for each washed IP (from Step 3 in Section D) and the 2% input sample (from Step 3 in Section C).

1. Add 150 µl of the 1X ChIP Elution Buffer to the 2% input sample tube and set aside at room temperature until Step 7.

2. Pellet Protein G Agarose Beads by brief 1 min centrifugation at 6000 rpm and remove supernatant.
3. Add 150 µl 1X ChIP Elution Buffer to each IP sample.
4. Elute chromatin from the antibody/Protein G beads for 30 minutes at 65°C with gentle vortexing (1,200 rpm). A thermomixer works best for this step. Alternatively, elutions can be performed at room temperature with rotation, but may not be as complete.
5. Pellet Protein G Agarose Beads by brief 1 min centrifugation at 6,000 rpm.
6. Carefully transfer eluted chromatin supernatant to a new tube.
7. To all tubes, including the 2% input sample from Step 1, reverse cross-links by adding 6 µl 5M NaCl and 2 µl Proteinase K, and incubate 2 h at 65°C. This incubation can be extended overnight.
8. Immediately proceed to Section F. Alternatively, samples can be stored at -20°C. However, to avoid formation of a precipitate, be sure to warm samples to room temperature before adding DNA Binding Reagent A (Section F, Step 1).

### F. DNA Purification Using Spin Columns:

**Before starting:**

- Add 12 ml of isopropanol to DNA Binding Reagent A and 24 ml of ethanol (96-100%) to DNA Wash Reagent B before use. These steps only have to be performed once prior to the first set of DNA purifications.
- Remove one DNA purification spin column and collection tube for each DNA sample from Section E.

1. Add 600 µl DNA Binding Reagent A to each DNA sample and vortex briefly.
  - 4 volumes of DNA Binding Reagent A should be used for every 1 volume of sample.
2. Transfer 375 µl of each sample from Step 1 to a DNA purification spin column in collection tube.
3. Centrifuge at 14,000 rpm in a microfuge for 30 sec.
4. Remove the spin column from the collection tube and discard the liquid. Replace spin column in the collection tube.
5. Transfer the remaining 375 µl of each sample from Step 1 to the spin column in collection tube. Repeat Steps 3 and 4.
6. Add 700 µl of DNA Wash Reagent B to the spin column in collection tube.
7. Centrifuge at 14,000 rpm in a microfuge for 30 sec.
8. Remove the spin column from the collection tube and discard the liquid. Replace spin column in the collection tube.
9. Centrifuge at 14,000 rpm in a microfuge for 30 sec.
10. Discard collection tube and liquid. Retain spin column.
11. Add 50 µl of DNA Elution Reagent C to each spin column and place into a clean 1.5 ml microfuge tube.
12. Centrifuge at 14,000 rpm in a microfuge for 30 sec to elute DNA.
13. Remove and discard DNA purification spin column. Eluate is now purified DNA. Samples can be stored at -20°C.

### G. Quantification of DNA by PCR:

**Recommendations:**

- Use Filter-tip pipette tips to minimize risk of contamination.
- The control primers included in the kit are specific for the human or mouse RPL30 gene and can be used for either standard PCR or quantitative real-time PCR. If the user is performing ChIPs from another species, it is recommended that the user design the appropriate specific primers to DNA and determine the optimal PCR conditions.
- A Hot-Start Taq polymerase is recommended to minimize the risk of non-specific PCR products.
- PCR primer selection is critical. Primers should be designed with close adherence to the following criteria:

## Chromatin Immunoprecipitation Protocol (cont.)

<b>Primer length:</b>	24 nucleotides
<b>Optimum Tm:</b>	60°C
<b>Optimum GC:</b>	50%
<b>Amplicon size:</b>	150 to 200 bp (for standard PCR) 80 to 160 bp (for real-time quantitative PCR)

### Standard PCR Method:

1. Label the appropriate number of 0.2 ml PCR tubes for the number of samples to be analyzed. These should include the 2% input sample, the positive control Histone H3 sample, the negative control Normal Rabbit IgG sample, and a tube with no DNA to control for DNA contamination.
2. Add 2 µl of the appropriate DNA sample to each tube.
3. Prepare a master reaction mix as described below, making sure to add enough reagent for two extra tubes to account for loss of volume. Add 18 µl of master mix to each reaction tube.

Reagent	Volume for 1 PCR Reaction (18 µl)
Nuclease-free H <sub>2</sub> O	12.5 µl
10X PCR Buffer	2.0 µl
4 mM dNTP Mix	1.0 µl
5 µM RPL30 Primers	2.0 µl
Taq DNA Polymerase	0.5 µl

4. Start the following PCR reaction program:

a.	Initial Denaturation	95°C	5 min
b.	Denature	95°C	30 sec
c.	Anneal	62°C	30 sec
d.	Extension	72°C	30 sec
e.	Repeat Steps b-d for a total of 34 cycles.		
f.	Final Extension	72°C	5 min

5. Remove 10 µl of each PCR product for analysis by 2% agarose gel or 10% polyacrylamide gel electrophoresis with a 100 bp DNA marker. The expected size of the PCR product is 161 bp for human RPL30 and 159 bp for mouse RPL30.

### Real-Time Quantitative PCR Method:

1. Label the appropriate number of PCR tubes or PCR plates compatible with the model of PCR machine to be used. PCR reactions should include the positive control Histone H3 sample, the negative control Normal Rabbit IgG sample, a tube with no DNA to control for contamination, and a serial dilution of the 2% input chromatin DNA (undiluted, 1:5, 1:25, 1:125) to create a standard curve and determine the efficiency of amplification.
2. Add 2 µl of the appropriate DNA sample to each tube or well of the PCR plate.
3. Prepare a master reaction mix as described below. Add enough reagents for two extra reactions to account for loss of volume. Add 18 µl of reaction mix to each PCR reaction tube or well.

Reagent	Volume for 1 PCR Reaction (18 µl)
Nuclease-free H <sub>2</sub> O	6 µl
5 µM RPL30 Primers	2 µl
SYBR-Green Reaction Mix	10 µl

4. Start the following PCR reaction program:
  - a. Initial Denaturation 95°C 3 min
  - b. Denature 95°C 15 sec
  - c. Anneal and Extension: 60°C 60 sec
  - d. Repeat steps b and c for a total of 40 cycles.
5. Analyze quantitative PCR results using the software provided with the real-time PCR machine. Alternatively, one can calculate the IP efficiency manually using the Percent Input Method and the equation shown below. With this method, signals obtained from each immunoprecipitation are expressed as a percent of the total input chromatin.

$$\text{Percent Input} = 2\% \times 2^{(C[T] \text{ 2\% Input Sample} - C[T] \text{ IP Sample})}$$

$C[T] = C_T = \text{Threshold cycle of PCR reaction}$

## APPENDIX A: Optimization of Chromatin Digestion

Optimal conditions for digestion of cross-linked DNA to 150-900 base pairs in length depend on cell type and cell concentration as well as the concentration of Micrococcal Nuclease. Below is a protocol to determine the optimal digestion conditions for a specific cell type and concentration of cells.

1. Prepare cross-linked nuclei from  $4 \times 10^7$  cells as described in Section A, Steps 1-6.
2. Transfer 200 µl of the nuclei preparation from Step 6 in Section A into 5 individual microcentrifuge tubes and place on ice.
3. Add 5 µl Micrococcal Nuclease stock to 20 µl of 1X Buffer B + DTT (1:5 dilution of enzyme).
4. To each of the 5 tubes in Step 2, add 0 µl, 2.5 µl, 5 µl, 7.5 µl or 10 µl of the diluted Micrococcal Nuclease, mix by inverting tube several times and incubate for 20 min at 37°C with frequent mixing.
5. Stop digest by adding 20 µl of 0.5 M EDTA and placing tubes on ice.
6. Pellet nuclei by centrifugation at 13,000 rpm in a microcentrifuge for 1 min at 4°C and remove supernatant.
7. Resuspend nuclear pellet in 200 µl of 1X ChIP buffer + PIC + PMSF. Incubate on ice for 10 min.
8. Sonicate lysate with several pulses to break nuclear membrane. Incubate samples 30 sec on wet ice between pulses.
  - Optimal conditions required for complete lysis of nuclei can be determined by observing nuclei under light microscope before and after sonication.
  - HeLa nuclei were completely lysed after 3 sets of 20-second pulses using a VirTis
9. Clarify lysates by centrifugation at 10,000 rpm in a microcentrifuge for 10 min at 4°C.
10. Transfer 50 µl of each of the sonicated lysates to new microfuge tubes.
11. To each 50 µl sample, add 100 µl nuclease-free water, 6 µl 5 M NaCl and 2 µl RNase A. Vortex to mix and incubate samples at 37°C for 30 min.
12. To each RNase A-digested sample, add 2 µl Proteinase K. Vortex to mix and incubate sample at 65°C for 2 hours.
13. Remove 20 µl of each sample and determine DNA fragment size by electrophoresis on a 1% agarose gel with a 1 kb DNA marker.
14. Observe which of the digestion conditions produces DNA in the desired range of 150-900 base pairs (1 to 5 nucleosomes, see Figure 1). The volume of diluted Micrococcal Nuclease that produces the desired size of DNA fragments using this optimization protocol is equivalent to the volume of Micrococcal Nuclease stock that will need to be added to  $4 \times 10^7$  cells to produce the desired size of DNA fragments. For example, if 5 µl of diluted Micrococcal Nuclease produces DNA fragments of 150-900 base pairs in this protocol, then 5 µl of stock Micrococcal Nuclease should be added to  $4 \times 10^7$  cells during the digestion of chromatin in Step 7 of Section A.
15. If results indicate that DNA is not in the desired size range, then repeat optimization protocol, adjusting the amount of Micrococcal Nuclease in each digest accordingly.

Virsonic 100 Ultrasonic Homogenizer/Sonicator set at setting 6 with a 1/8-inch probe.

- Alternatively, nuclei can be lysed by homogenizing the lysate 20 times in a Dounce homogenizer; however, lysis may not be as complete.

## APPENDIX B: Troubleshooting Guide

Problem	Possible Causes	Recommendation
1. Concentration of the digested chromatin is too low.	Not enough cells added to the chromatin digestion or nuclei were not completely lysed after digestion.	If DNA concentration of the chromatin preparation is close to 50 µg/ml, add additional chromatin to each IP to give at least 5 µg/IP and continue with protocol.  Count a separate plate of cells before cross-linking to determine an accurate cell number and/or visualize nuclei under microscope before and after sonication to confirm complete lysis of nuclei.
2. Chromatin is under-digested and fragments are too large (greater than 900 bp).	Cells may have been over cross-linked. Cross-linking for longer than 10 min may inhibit digestion of chromatin.  Too many cells or not enough Micrococcal Nuclease was added to the chromatin digestion.	Perform a time course at a fixed formaldehyde concentration. Shorten the time of cross-linking to 10 min or less.  Count a separate plate of cells before cross-linking to determine accurate cell number and see Appendix A for optimization of chromatin digestion.
3. Chromatin is over-digested and fragments are too small (exclusively 150 bp mono-nucleosome length). Complete digestion of chromatin to mono-nucleosome length DNA may diminish signal during PCR quantification, especially for amplicons greater than 150 bp in length.	Not enough cells or too much Micrococcal Nuclease added to the chromatin digestion.	Count a separate plate of cells before cross-linking to determine accurate cell number and see Appendix A for optimization of chromatin digestion.
4. No product or very little product in the input PCR reactions.	Not enough DNA added to the PCR reaction or conditions are not optimal.  PCR amplified region may span nucleosome-free region.  Not enough chromatin added to the IP or chromatin is over-digested.	Add more DNA to the PCR reaction or increase the number of amplification cycles.  Optimize the PCR conditions for experimental primer set using purified DNA from cross-linked and digested chromatin. Design a different primer set and decrease length of amplicon to less than 150 bp (see primer design recommendations in Section H).  For optimal ChIP results add 5–10 µg chromatin per IP. See recommendations for problems 1 and 3 above.
5. No product in the positive control Histone H3-IP RPL30 PCR reaction.	Not enough chromatin or antibody added to the IP reaction or IP incubation time is too short.  Incomplete elution of chromatin from Protein G beads.	Be sure to add 5–10 µg of chromatin and 10 µl of antibody to each IP reaction and incubate with antibody over-night and an additional 2 h after adding Protein G beads.  Elution of chromatin from Protein G beads is optimal at 65°C with frequent mixing to keep beads suspended in solution.
6. Quantity of product in the negative control Rabbit IgG-IP and positive control Histone H3-IP PCR reactions is equivalent.	Too much or not enough chromatin added to the IP reaction. Alternatively, too much antibody added to the IP reaction.  Too much DNA added to the PCR reaction or too many cycles of amplification.	Add no more than 15 µg of chromatin and 10 µl of histone H3 antibody to each IP reaction. Reduce the amount of normal rabbit IgG to 1 µl per IP.  Add less DNA to the PCR reaction or decrease the number of PCR cycles. It is very important that the PCR products are analyzed within the linear amplification phase of PCR. Otherwise, the differences in quantities of starting DNA can not be accurately measured.
7. No product in the Experimental Antibody-IP PCR reaction.	Not enough DNA added to the PCR reaction.  Not enough antibody added to the IP reaction.  Antibody does not work for IP.	Add more DNA to the PCR reaction or increase the number of amplification cycles.  Typically a range of 1 to 5 µg of antibody are added to the IP reaction; however, the exact amount depends greatly on the individual antibody. Increase the amount of antibody added to the IP.  Find an alternate antibody source.