

qPCR 8-STRIP COMPATIBILITY CHART

	MB-q100	MB-q200	Note/Holder NO.
Standard Well (0.2ml, height 21mm.)		☀	
Low Profile Well (0.1ml, height 16mm.)	☀		
Applied Biosystems			
7000		☀	
7300		☀	Precision Plate Holder for 0.2ml tubes and strips(4367033)
7500		☀	Precision Plate Holder for 0.2ml tubes and strips(4367033)
StepOne	☀		MicroAmp™ Fast 48-Well Tray (4375282)
StepOne Plus	☀		MicroAmp™ 96-Well Tray for VeriFlex™ Blocks (4379983)
7500Fast	☀		Precision Plate Holder for 0.1ml tubes and strips(4403809)
7900HT Fast (with tray) (0.1ml type)			Use Plate Only
7900HT (with tray) (0.2ml type)			Use Plate Only
ViiA7 Fast(0.1ml model)	☀		ViiA™ 7 Fast 96-Well Tube Adaptor(4462078) & Retainer Tray for Fast 96-well (0.1mL) tubes and strips(4379983)
ViiA7 (0.2ml model)		☀	ViiA™ 7 96-Well Tube Adaptor(4459845) & Retainer Tray for Standard 96-well (0.2mL) tubes and strips(4381850)
"QuantStudio Fast 96-well Block(0.1ml model)"	☀		MicroAmp® 96-Well Tray (Black) (for 0.1ml)(4379983)
"QuantStudio 96-well Block(0.2ml model)"		☀	MicroAmp® 96-Well Tray/Retainer Set (Blue) (for 0.2ml) (4381850)
Roche			
LightCycler 480			
LightCycler 96	☀		
LightCycler Nano	☀		32-well
Bio-Rad			
iCycler iQ		☀	
iQ5		☀	
MyiQ		☀	
Chromo4		☀	
Opticon I	☀		
Opticon II	☀		
mini Opticon	☀		48-Well
CFX96	☀		
CFX384			
Stratagene			
Mx4000/3000P/3005P		☀	
Eppendorf			
MasterCycle ep Realplex		☀	
Thermo Scientific			
PikoReal 24 Real-Time PCR	☀		24-Well
Biometra			
Topical Thermocycler		☀	